
PH-31

MRI QA Phantom

Instruction
Manual

Contents

- Before beginning
 - Product contents and warning P.1
 - How to use the phantom P.2-P.3
- Preparation
 - Handling the phantom P.4
 - About MRI contrast solution
 - How to fill with MRI solution P.5-P.6
 - How to fill with MRI solution - Warning P.7
- Resources P.8-P.9

Product contents

Before your first use, ensure that your product includes all components listed below.

a	Phantom unit A	1	e	MRI contrast solution	5
b	Phantom unit B	1		*5 levels of concentration, 50mL per bottle	
c	MRI contrast solution (10mM 10L)	1	f	Petroleum jelly	1
d	Spout for MRI contrast solution	1	g	Screwdriver	1
			h	Extra screws	4
				Instruction manual	

Warning

- Please handle with care. The phantom is made of hard resinous material; the phantom can be damaged if dropped or hit by a hard object.

- Please clean the phantom with water or pH neutral detergent. Please do not use organic detergent such as thinner.

- Please do not store in high temperatures, high humidity or in direct sunlight. It may lead to deformation or damages.

- Please do not write with a pen. Writing directly on the phantom with felt-tip pen or permanent marker will not wash off.


How to use the phantom

Phantom unit A

Liquid capacity 3180mL Acrylic materials, Outer diameter 220mm


Structural design


How to use the phantom

Phantom unit B

Liquid capacity 3300mL Acrylic materials, Outer diameter 220mm

Structural design

1 Handling the phantom

The recommended temperature for use of this phantom is 4 degrees above or below 22 degrees Celcius according to the NEMA and AAPM standards. For the phantom to be at the same temperature as the scan room, place the phantom in the room for 24 hours before scanning.

MRI signal materials are sensitive to temperature change. For every 1 degree Celcius of increase, T1 will increase by approximately 1% and T2 will decrease by approximately 1%.

Additionally, the phantom is composed of plastic material and is prone to damage from excessive force or reactive solvent.

2 About MRI contrast solution

The MRI is one of the most advanced medical imaging diagnosis method. However, standard uses and materials for precise product management and maintenance have not yet been defined.

The MRI contrast liquid included in this product has been developed with a focus in safety and has been featured in the Report of AAPM*. The primary ingredient is nickel dichloride (NiCl_2).

The table on page 15 displays the values of T1 and T2 of the magnetic resonance of the nickel dichloride. Figure 1 on page 16 is a graph of the data on page 15.

Kyoto Kagaku's phantom includes contrast solution of 5 different levels of concentration to compare with the 10mM solution as a standardization product.

*Quality assurance methods and phantoms for magnetic resonance imaging


: Report of AAPM nuclear magnetic resonance Task Group No.1

3 How to fill with MRI solution

1. Remove the cap on the MRI contrast liquid container, and attach the spout on the opening. Prepare a container (ie. beaker). Open the spout and carefully tilt the box to pour the solution into the container.

2. Loosen and remove the screws of the phantom lid with the screw driver. Do not remove the white plugs at this time. Please take care not to damage the screws when unscrewing.


3 How to fill with MRI solution

3. Carefully pour the nickel chloride solution into the phantom, filling up to 50% to 80% of the body.

At this point, apply a small amount of jelly to the rubber portion between the lid and the main body of the phantom.


4. Once the lid is attached back on the main body of the phantom, tilt the phantom to remove any air bubbles in the solution.

Removed plug


5. Remove the white plugs and pour in the contrast solution through the opening to completely fill the phantom.

Caution: Please use the "V" mark to align the lid and the main body of the phantom.

3 How to fill with MRI solution

6. Insert one plug. Tilt the phantom and remove any remaining air bubbles in the solution.

7. Insert the second plug. Be sure to not add any air bubbles.

4 Warning


T1 and T2 values are heavily dependent on the temperature change; please take extreme care in maintaining a consistent temperature when the phantom is in use. Leave the phantom with the contrast solution in room temperature for over 24 hours to retain a consistent temperature for the product. Raw materials used to manufacture the provided nickel dichloride solution is manufactured in 23 degrees Celcius with relative humidity of 40%.

Resources

NiCl ₂ density (mM)	1	5	6	7	8
T1 value (msec)	1117	329	264	229	203
Standard deviation	23.53	5.54	10.11	4.5	3.02
T2 value (msec)	1044	291	246	213	189
Standard deviation	6.96	0.83	0.56	0.61	0.66
NiCl ₂ density (mM)	9	10	15	20	25
T1 value (msec)	184	174	116	89	69
Standard deviation	4.28	2.11	1.27	1.03	0.71
T2 value (msec)	169	152	102	79	64
Standard deviation	0.51	0.65	0.5	0.11	0.4

T1 and T2 values of the magnetic resonance of varying nickel dichloride densities
*Data from room temperature 24 degree Celcius and magnetic field strength of 1.5 Tesla

Resources

T1 and T2 values of varying nickel dichloride densities

URL : <http://www.kyotokagaku.com>
e-mail : rw-kyoto@kyotokagaku.co.jp

Worldwide Inquiries & Ordering

Kyotokagaku Head Office and Factories:

15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN
Tel: +81-75-605-2510 Fax: +81-75-605-2519

Kyotokagaku USA Office: USA and Canada sales and services

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA
Tel: 1-310-325-8860 Fax: 1-310-325-8867

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