Scrotal Ultrasound Phantom

Master techniques for a safe and painless examination


This phantom helps accurate diagnosis of scrotum abnormalities through hands-on training

Product Supervision

Charlotte Henningsen, MS, RT(R), RDMS, RVT, ESDMS, FAIUM Chair & Professor - Sonography Department Adventist University of Health Sciences

FEATURES

- 1 | Two types of interchangeable scrotums: normal and pathology
- 2 | Excellent ultrasound image quality



SKILLS


| Patient positioning and preparation for examination


- Scrotal ultrasound screening
- | Visualization of testicular cancer



Normal phantom

- scrotum, testicle, epididymis and penis









Pathological phantom

- -scrotum, testicle, epididymis and penis
- -testicular cancer (each one in left and right testis, 10mm dia.)







DESCRIPTIONS

SPECIFICATIONS

Size: 34x33x24cm/13.4x13x9.4 inch Weight: 3.6kg/7.9lbs

MATERIALS

- 1 pelvic body
- normal scrotum
- 1 pathological scrotum
- storage case
- 1 instruction manual

IREPLACEMENT PARTS

41911-010 1 normal scrotum unit 41911-020 1 pathological scrotum unit

* Ultrasound scanner is not included.





