

# Mobile Holder for Image Receptors

Cat# MCH-33B

1.5" Deep Channel Accepts Grids and Caps and will accommodate most receptor sizes like: Agfa, AllPro, Carestream, Canon, Dell, Fuji, GE, Kodak, Konica, Quantum, Medlink, Rayence, Samsung, Shimadzu, Toshiba, Varian, and Vieworks

Black arm is 12"

Front of cassette is 18" from pole

60" Chromed 2" x 2" Steel


Posi-Lock height adjustable arm with spring-loaded pin.



## Tilt & Rotation Abilities:


- \* Straight Vertical +/- 90 degrees
- \* Incline Forward / Backward
- \* InclHorizontal rotation +/- 180 degrees
- \* Full Rotate or Oblique Tilt
- \* Face Up or Down
- \* 360 degrees - Anywhere in Between

Vertical Range from bottom of image receptor: 1" off Floor to 54"



Sturdy 24" x 24" Broad H-base prevents tipping.

Heavy Duty locking 3" casters for easy manipulation.



- Platinum baked enamel finish. Easy to clean base is perfect for OR settings.
- Bilateral clamp safely secures imaging receptors in place. Tilt & Rotate Range of Motion allows for positioning in lateral, horizontal, vertical and oblique positions.
- Channel is 1.5" deep to accept grids and caps
- Will accept receptor sizes up to 20" tall. No width limitation.
- Holds receptors up to 12 lbs.
- Manufactured in Nashville, Tennessee USA
- 14 Day Lead Time
- One Year Warranty
- \* Ships in 3 cartons: Pole 12lbs, base 23 lbs, & arm 19 lbs

