



Mobile Holder for Image Receptors

Cat# MCH-37B



Tilt & Rotation Abilities:

- * Straight Vertical +/- 90 degrees
- * Incline Forward / Backward
- * Horizontal rotation +/- 180 degrees
- * Full Rotate or Oblique Tilt
- * Face Up
- * Face Down
- * 360 degrees - Anywhere in Between



Image Receptor Obliqued at Angle



Image Receptor Face Up / Down



Image Receptor Tilted Backwards

- Radiographs can be taken from any angle. Tray can be positioned face up, face down, and anywhere in between to maximize proper receptor positioning. Perfect for ER, OR, and x-ray procedures such as myelograms, decubitus studies, and three way acute abdomen exams.
- Platinum baked enamel finish base that is easy to clean base is perfect for OR
- Will hold receptors up to 12 lbs.
- Will accept flat image receptors up to 20" tall and no width limitations.
- Manufactured in Nashville, Tennessee USA
- One Year Warranty
- 14 Day Lead Time
- * Ships in 3 cartons: Pole 12lbs, base 23 lbs , & arm 19 lbs.

Heavy Duty locking 3" casters for easy manipulation.

