


# Mobile Holder for Image Receptors

Cat# MCH-37C



1.5" Deep Channel

Telescoping Extension Arm Extends 11" to 23" from

Vertical Range from bottom of image receptor: Floor to 54" Allowing ample room for imaging at several height levels.



Sturdy 24" x 24" Broad H-base prevents tipping.

Heavy Duty locking 3" casters for easy manipulation.

1.5" Deep Channel Accepts Grids and Caps and will accommodate most receptor sizes like: Agfa, AllPro, Carestream, Canon, Dell, Fuji, GE, Kodak, Konica, Quantum, Medlink, Rayence, Samsung, Shimadzu, Toshiba, Varian, and Vieworks



- \* Bilateral Clamps securely holds receptors in place.
- \* Holds up to 20.5" wide imaging receptors. No Height Limitations.
- \* Will hold receptors up to 12 lbs.
- \* Easy to clean in a "OR" or "Sterile" setting.
- \* Manufactured in Nashville, Tennessee USA
- \* 14 Day Lead Time
- \* One Year Warranty
- \* Ships in 3 cartons: Pole 12lbs, base 23 lbs , & arm 19 lbs.



- Tilt & Rotation Abilities:**
- \* Straight Vertical
  - \* Incline Forward
  - \* Incline Backward
  - \* Rotate or Oblique Tilt
  - \* Face Up
  - \* Face Down
  - \* Anywhere in Between