

Mobile Holder for Image Receptor

Cat# MCH-48B



- Tall Format Lateral (complete mobile unit)
- Holds 14" x 36" or 51" image receptors for AP/PA images. Fixed position only. This is non-articulating).
- This will not rotate to accommodate Cross table Lateral studies.
- Mobile unit includes a vertical range from bottom of image receptor: floor to 54".
- Posi-lock height adjustable arm with spring loaded pin.
- Features 60" chromed steel post. 24" H-Base with locking casters is broad to prevent tipping.
- Platinum baked enamel finish. Easy to clean base is perfect for OR settings.
- Bilateral clamp safely secures imaging receptors in place.
- Channel is 1.5" deep to accept grids and caps.
- * Manufactured in Nashville, Tennessee USA
- * 14 Day Lead Time
- * Ships in three cartons: Pole, base, and arm
- * For Cross-Table Lateral studies we recommend to order MCH-45B

