

# MR Height Adjustable Trolley (Gen 2)

## The Adjustable and MRI Safe Patient Transportation Solution



(7 Tesla)

Imaging Solutions' new MRI-safe, height adjustable trolley provides transportation to and from the MRI room with easy manoeuvrability. Constructed from non-ferrous materials, this patient trolley is usable in environments up to 7.0T.

The new patient trolley features a fully adjustable height and backrest, as well as patented Trendelenburg tilt technology which makes for quick operation during emergency situations. The trolley is durable, easy to use, and ideal for any hospital.



### Height Adjustable

Easy height adjustment with dual pedal pneumatic operation supporting major MRI vendors.



### Trendelenburg Tilt

Imaging Solutions patented Trendelenburg tilt technology for quick operation during emergency situations.



### MRI Compatible

Constructed from non-ferrous materials, usable in environments up to 7.0 T.



### Australian Made

Designed and manufactured in Australia offering complete peace of mind. Local expertise and components when you need them.



### Adjustable Backrest

Discreet friction lock allows easy adjustment to any angle.



### Direction Locking / Central Braking

Improved mobility through innovative directional locking. Enhanced safety with full castor locking with a single motion.

your single source supplier



### Side Rails

Unobtrusive side rail design disappears when not required allowing you to get closer to the bed or navigate those tight corridors.



### TENTE® Mobility

Using medical-grade castors from market leader TENTE. Featuring non-marking tread and precision ball bearings, you can be assured of superior mobility performance.



### Corner Buffering Protection

Protect your trolley and your facility from unintended collisions with soft plastic corner caps.



### AEGIS® Mattress

Our mattresses are manufactured using premium grade sterilised foam and coated with our AEGIS medical grade fabric allowing for sterilisation, durability and patient comfort.



### Durable Finish

Sandblasted and powder coated finish ensuring your trolley stays newer for longer.



### Included Warranty

Unbeatable five year structural warranty on all parts and expert workmanship will provide you peace of mind for years to come.



### Unique Design

Simple, innovative and efficient design with fewer potential points of failure.

## Specifications

<b>Overall Dimensions</b>	2042 x 747 mm (w/ side rails up) 2042 x 657 mm (w/ side rails down)
<b>Patient Surface</b>	1800 x 550 mm (backrest 750 mm)
<b>Height Range</b>	530 to 900 mm (without mattress)
<b>Weight Capacity</b>	250kg
<b>MR Compatibility</b>	Up to 7.0 T
<b>Side Rail Height</b>	135 mm (above mattress)
<b>Trendelenburg Tilt</b>	15 degrees
<b>Mattress Thickness</b>	50 mm

## Optional Accessories

- IV Pole Mounts
- Oxygen Cylinder Holder
- Mattresses (customised or replacement)
- Optional Frame Colours Available

## Product Code

Height Adjustable Trolley	ISM-HT-58123
Fixed Height Trolley	ISM-HT-58122



**Australia**  
Phone: 1300 132 100  
Fax: 1300 721 850

**New Zealand**  
Phone: 0800 723 776  
Fax: +61 7 3209 9812

**Service Hotlines**  
Orders: 1300 132 100  
Service: 1800 300 100

**Internet**  
www.imgsol.com  
info@imgsol.com