

MR Fixed Height Trolley (Gen 2.1)

The MRI Safe Patient Transportation Solution



(7 Tesla)

Imaging Solutions' new MRI-safe, fixed height trolley provides transportation to and from the MRI room with easy manoeuvrability. Constructed from non-ferrous materials, this patient trolley is usable in environments up to 7.0T.

The trolley is durable, easy to use, and ideal for any hospital.



MRI Compatible

Constructed from non-ferrous materials, usable in environments up to 7.0 T.



Australian Made

Designed and manufactured in Australia offering complete peace of mind. Local expertise and components when you need them.



Adjustable Backrest

Discreet friction lock allows easy adjustment to any angle.



Direction Locking / Central Braking

Improved mobility through innovative directional locking. Enhanced safety with full castor locking with a single motion.



Side Rails

Unobtrusive side rail design disappears when not required allowing you to get closer to the bed or navigate those tight corridors.




AEGIS® Mattress

Our mattresses are manufactured using premium grade sterilised foam and coated with our AEGIS medical grade fabric allowing for sterilisation, durability and patient comfort.

your single source supplier





Corner Buffering Protection

Protect your trolley and your facility from unintended collisions with soft plastic corner caps.



TENTE® Mobility

Using medical-grade castors from market leader TENTE. Featuring non-marking tread and precision ball bearings, you can be assured of superior mobility performance.



Durable Finish

Sandblasted and powder coated finish ensuring your trolley stays newer for longer.



Included Warranty

Unbeatable five year structural warranty on all parts and expert workmanship will provide you peace of mind for years to come.



Unique Design

Simple, innovative and efficient design with fewer potential points of failure.

Specifications

Overall Dimensions	2042 x 747 mm (w/ side rails up) 2042 x 657 mm (w/ side rails down)
Patient Surface	1800 x 550 mm (backrest 750 mm)
Weight Capacity	250kg
MR Compatibility	Up to 7.0 T
Side Rail Height	135 mm (above mattress)
Mattress Thickness	50 mm

Optional Accessories

- IV Pole Mounts
- Oxygen Cylinder Holder
- Mattresses (customised or replacement)
- Optional Frame Colours Available

Product Code

Height Adjustable Trolley	ISM-HT-58123
Fixed Height Trolley	ISM-HT-58122



Australia
Phone: 1300 132 100
Fax: 1300 721 850

New Zealand
Phone: 0800 723 776
Fax: +61 7 3209 9812

Service Hotlines
Orders: 1300 132 100
Service: 1800 300 100

Internet
www.imgsol.com
info@imgsol.com