
Accessory Product Guide

Fine-tune carts and cabinets with the accessories you need, located where you need them.

Attached Accessories page 2

Accessories attached to the exteriors of carts keep essential tools within easy reach, saving time and providing easy access to essential supplies and equipment. Sloped tops complete the attached accessory offering.

Procedure-specific Accessories page 9

Accessories target specific procedures and departments such as endoscopy and cardiac cath labs.

Interior Storage page 11

A modular system of shelves, dividers, trays, and baskets lets you arrange the interior configuration of your cabinets and carts for specific types and sizes of supplies and easily reconfigure as supplies and needs change.

Locks page 19

Select the type of security you prefer with keyed and keyless options, including proximity locks.

Casework Finishing Accessories page 20

Add countertops and finishing touches to the installation of Evolve casework.

Side-Mounted Accessory Rail (SAMR)

■ ROVER ■ TEMPO ■ ROAM



Accessory Rail attaches to the side of a Tempo cart, Roam cart, or Rover workstation to hold a variety of accessories. Attachment hardware included.

Dimensions

- 13.50"w x 1.50"h (for Rover workstation)
- 15.75"w x 1.50"h (for Tempo cart)
- 23.50"w x 1.50"h (for Roam cart)

Top-Mounted Accessory Rail (SUAR1, SUAR2, SNUAR1, SNUAR2)

■ TEMPO



Accessory rail attaches to the top of a Tempo cart to hold a variety of accessories. Single rail has 1 rail with posts; double rail has 2 rails with posts. Maximum rail height is 24" above the cart surface. Attachment hardware included.

- Specify SUAR1 for Accessory Rail with 1 rail
- Specify SUAR2 for Accessory Rail with 2 rails
- Specify SNUAR1 for Accessory Rail with 1 rail for Narrow Tempo cart
- Specify SNUAR2 for Accessory Rail with 2 rails for Narrow Tempo cart

Dimensions

- 27.50"w x 4.25"d (accessory rail only)

Breakaway Lock Bar (STLB)

■ TEMPO



Attaches to a Tempo cart to provide additional security of cart contents. The lock bar is sized to the length of the Tempo cart. Included with the Tempo Crash Cart.

Chart Holder (SCH)

■ ROVER ■ TEMPO ■ ROAM



Wire Chart Holder attaches to the side rail on a Tempo cart, Roam cart, or Rover workstation. Holds up to 5 pounds. Attachment hardware included.

Dimensions

- 13.25"w x 3.25"d x 7.75"h

Cord Wrap (SCW)

■ ROVER ■ TEMPO



Attaches to the back of a Tempo cart or Rover workstation. Attachment hardware included. Does not include power strip (see page 5).

Dimensions

- 12"w x 2.5"d x .75"h

Defibrillator Shelf (SDS)

■ TEMPO



Attaches to the top of a Tempo cart and holds up to 25 pounds. Cannot be used with Narrow Tempo cart. Comes with attachment hardware. Included with Tempo Crash Cart.

Dimensions

- 17"w x 13.75"d
- 14" clearance between top of cart and bottom of shelf

Embolectomy Storage Rack (SESR)

■ ROAM



Metal storage rack attaches to a fixed side rail on a Roam cart to organize catheters. Holds up to 25 pounds. Cannot be used on cart with an adjustable side-mounted rail. 5 fixed dividers and attachment hardware included.

Dimensions

- 18"w x 3.5"d x 30"h
- 2.93" space between dividers

Flip-Up Work Surface (SFWS)

■ ROVER ■ TEMPO ■ ROAM



Attaches to a fixed side rail on a Tempo cart, Roam cart, or Rover workstation. Work surface is made of ABS thermoplastic polymer and has seamless construction for cleaning ease. Holds up to 30 pounds. Cannot be used with Narrow Tempo cart or cart with an adjustable side-mounted rail. Attachment hardware included.

Dimensions

- 17"w x 13.75"d

IV Pole (SIVP)

■ TEMPO



Adjustable IV pole attaches to the top of a Tempo cart to hold and secure 2 IV bags. Adjusts from 34" to 56" and holds up to 10 pounds. Included with Tempo IV Prep Cart and Crash Cart.

Glove Box Holder (S1GBH, S3GBH)

■ ROVER ■ TEMPO ■ ROAM



Attaches to the side rail on a Tempo cart, Roam cart, or Rover workstation. Holds up to 5 pounds. Wire holder comes in 2 sizes and includes attachment hardware. Included with Tempo IV Prep Cart.

- Specify S1GBH for single box holder
- Specify S3GBH for triple glove box holder

Dimensions

- 11.5"w x 4"d x 4.25"h (single)
- 11.5"w x 4"d x 14.25"h (triple)

LED Task Light (STTL)

■ TEMPO



Adjustable LED task light attaches to the top of a Tempo cart. The light has a 26"-long flexible neck and includes a bulb with 2- and 5-watt settings and attachment hardware.

Monitor Arm (SFPMA)

■ TEMPO



Attaches to the top of a Tempo cart to hold monitors up to 20 pounds and 27" in width. Offers a height range of up to 25.75" and full tilt motion. Cannot be used with Narrow Tempo cart. Attachment hardware included.

Oxygen Tank Holder (SOTH)

■ ROVER ■ TEMPO ■ ROAM



Attaches to a fixed side rail on a Tempo cart, Roam cart, or Rover workstation to hold tanks up to 4" in diameter. Cannot be used on cart with an adjustable side-mounted rail. Holds up to 12 pounds and accommodates A, B, C, D, and E type tanks. Comes with attachment hardware. Included with Tempo Crash Cart.

Dimensions

- 6"w x 5"d x 19.50"h

Power Strip (SPS)

■ ROVER ■ TEMPO



For use with a Tempo cart or Rover workstation. 120-volt medical-grade power strip has 6 outlets, on/off switch, and 15' cord. Does not include cord wrap (see page 3).

Dimensions

- 13.75"w x 1.75"d x 2.5"h

Pull-Out Keyboard Tray (SKBT)

■ TEMPO



Fits into the top of a Tempo cart's interior in place of a 9"-high drawer to hold standard-sized keyboard. Holds up to 10 pounds and has a 6" height-adjustment range. Cannot be used in Narrow Tempo cart. Mouse pad and attachment hardware included.

Dimensions

- 21"w x 10.5"d

Push Handle (SAPH)

■ ROVER ■ TEMPO



Stainless steel handle attaches to a fixed side rail on a Tempo cart or the back of a Rover workstation. Cannot be used on cart with an adjustable side-mounted rail. Attachment hardware included.

Dimensions

- 16.5"w x 2.25"d x .75"h, 1.5" interior depth (for Roam cart)
- 15.25"w x 2"d x .75"h, 1.25" interior depth (for Tempo cart)
- 14"w x 2"d x .75"h, 1.25" interior depth (for Rover workstation)

Push Handle on Tempo cart shown

Resuscitation Board (SRB)

■ TEMPO



Attaches to the back of a Tempo procedure cart. Cannot be used with Narrow Tempo cart. Attachment hardware included. Included with Tempo Crash Cart.

Dimensions

- 15"w x 20"h

Sharps Container (SSC)

■ ROVER ■ TEMPO ■ ROAM



Locking 5-quart Sharps Container attaches to a fixed side rail on a Tempo cart, Roam cart, or Rover workstation. Holds up to 25 pounds. Cannot be used on cart with an adjustable side-mounted rail. Comes with attachment hardware. Included with Tempo Anesthesia Cart and IV Prep Cart.

Dimensions

- 11.25"w x 6.25"d x 11.75"h

Side Storage Bin (SSSB)

■ ROVER ■ TEMPO ■ ROAM



Metal Storage Bin attaches to the side rail on a Tempo cart, Roam cart, or Rover workstation. Holds up to 5 pounds.

Dimensions

- 9.5"w x 5.5"d x 6"h

Sloped Top for Evolve Cabinet (THST, TRST, UHST, SHST, SRST)

■ EVOLVE



Adds a sloped top to an Evolve cabinet with hinged or roll-top door.
Attached hardware included.

- 9.75" h sloped top, for 27.25" d cabinet:
 - Specify 9THST for cabinet with hinged doors
 - Specify 9TRST for cabinet with roll-top door
- 7" h sloped top, for 19.25" d cabinet:
 - Specify 7THST for cabinet with hinged doors
 - Specify 7TRST for cabinet with roll-top door
- 4.75" h sloped top, for 13" d cabinet:
 - Specify 4UHST for cabinet with hinged doors
- 3" h vented sloped top, for scope cabinet:
 - Specify 3SHST for cabinet with hinged doors
 - Specify 3SRST for cabinet with roll-top door

Sloped Top for Roam Cart (ST)

■ ROAM



Adds a powder-coated 5"-high steel sloped top to a Roam cart with hinged door.
Attachment hardware included.

- Integrated sloped tops come with Roam carts with roll-top doors

Surface Rail (SWAR, SNWAR)

■ TEMPO



Aluminum surface rail attaches to the side rails on a Tempo cart to keep supplies and equipment from falling off surface. Attachment hardware included.

- Specify SWAR for use with Tempo Cart
- Specify SNWAR for use with Narrow Tempo Cart

Dimensions

- 27.25" w x 20.50" d x 5.75" h
- 4" clearance between top of cart and bottom of rail

Tilt Bin Module

■ TEMPO



Modules come in 1-, 4-, 5-, or 6-bin configurations and attach to the upper rail on a Tempo cart. Bin module is enclosed with clear plastic fronts. 4- and 5-bin module included with Tempo Anesthesia Cart; 4-bin module included with Tempo IV Prep Cart.

- Specify STB1 for 1-compartment Tilt Bin, for use with Narrow Tempo cart only
- Specify STB4 for 4-compartment Tilt Bin
- Specify STB5 for 5-compartment Tilt Bin
- Specify STB6 for 6-compartment Tilt Bin

Notes

- 4-, 5- and 6-bin modules cannot be used with Narrow Tempo cart
- Accessory Rail sold separately (SUAR1, SUAR2, SNUAR1, SNUAR2)

Dimensions

- STB1 – 5.5" x 5"w x 6.5"h
- STB4 – 6.5"d x 23.5"w x 8"h
- STB5 – 5.5"d x 23.5"w x 6.5"h
- STB6 – 3.5"d x 23.5"w x 4.5"h

Waste Basket (SWB)

■ ROVER ■ TEMPO ■ ROAM



Plastic Waste Basket attaches to a fixed side rail on a Tempo cart, Roam cart, or Rover workstation. Holds up to 15 pounds. Cannot be used on cart with an adjustable side-mounted rail. Attachment hardware included. Included with Tempo Anesthesia Cart.

Dimensions

- 11.5"w x 9"d x 13.25"h

Wire Basket (SSMWB)

■ ROVER ■ TEMPO ■ ROAM



White wire basket attaches to the side rail on a Tempo cart, Roam cart, or Rover workstation and holds up to 5 pounds. Attachment hardware included.

Dimensions

- 12"w x 6"d x 6"h

Cath Manager® (S7CM, S9CM)

■ ROAM ■ EVOLVE



Attaches to a Metal Shelf for Cath Managers in a Roam cart or Evolve cabinet with FlexCell interior to manage catheters. Each inline hook holds 10-15 catheter packages and can be adjusted horizontally along shelf. Uses 3 FlexCells. Includes product labeling.

- Specify S7CM for 7-hook cath manager; attaches to 1 or 2 compartment, 19.25"d shelf
- Specify S9CM for 9-hook cath manager; attaches to 1 or 2 compartment, 27.25"d shelf

Cath Manager Module (S19M1CM4, S19M2CM7, S27M1CM3, S27M2CM7)

■ ROAM ■ EVOLVE



Module includes a metal shelf and cath managers for use in a Roam cart or Evolve cabinet. Uses 3 FlexCells.

- Specify S19M1CM4 for 4 7-hook cath manager module; attaches to 1 compartment, 19.25"d shelf
- Specify S19M2CM7 for 7 7-hook cath manager module; attaches to 2 compartment, 19.25"d shelf
- Specify S27M1CM3 for 3 9-hook cath manager module; attaches to 1 compartment, 27.25"d shelf
- Specify S27M2CM7 for 7 9-hook cath manager module; attaches to 2 compartment, 27.25"d shelf

Fan Ventilation Unit (SFVU)

■ ROAM ■ EVOLVE



For use in all Evolve scope cabinets and a Roam scope cart. Includes HEPA filter and 10 AMP 125VAC 10' power cord.

Notes

- Order replacement filters separately (SHEPFIL for 1 filter or SHEPFIL3 for 3 filters)

Dimensions

- 7"w x 5"d x 5"h

Scope Manager® (SRSH1, SRSH2)

■ ROAM ■ EVOLVE



For use in all Evolve scope cabinets and a Roam scope cart. Includes metal bracket.

- Specify SRSH1 for single Scope Manager
- Specify SRSH2 for double Scope Manager

Side-Mounted Scope Cabinet (SSMSC)

■ TEMPO



Attaches to a fixed side rail on a Tempo cart to hold scopes up to 39" in length. The cabinet has a vented back, hinged acrylic door with keyed lock, and includes 3 scope managers with cord retainers. Cannot be used with Narrow Tempo cart or cart with an adjustable side-mounted rail. Attachment hardware included.

Dimensions

- 14.5"w x 12"d x 48"h

Suture Module (SPSM)

■ ROAM ■ EVOLVE



Includes 6 pull-out shelves to store and access suture boxes and other small supplies in a Roam cart or Evolve cabinet. Each shelf holds 8 2.25" deep suture boxes for a total of 48 boxes. Uses 11 FlexCells.

Notes


- Order suture shelf dividers separately (SSSHDI)
- Order additional shelves separately (SSSH)

Dimensions

- 14"w x 23.5"d x 26.5"h

4"H Basket, Clear (SB4)

■ ROAM ■ EVOLVE



For use in a Roam cart or Evolve cabinet with FlexCell interior. Perforated bottoms and sides keep dust out of baskets. Can be easily removed to clean cart or cabinet interior. Safety stops prevent basket from coming out of cart or cabinet. Baskets can be organized using adjustable dividers; order separately. Holds up to 40 pounds. Uses 2 FlexCells.

Notes

- Specify SBD4 for short adjustable divider
- Cell use with angled tilt: 5 FlexCells on 19.25"d cabinet and 7 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 23.5"w x 15.75"d x 4"h
- Interior: 21.5"w x 13.75"d x 3.75"h

■ ROAM ■ EVOLVE

8"H Basket, Clear (SB8)



For use in a Roam cart or Evolve cabinet with FlexCell interior. Perforated bottoms and sides keep dust out of baskets. Can be easily removed to clean cart or cabinet interior. Safety stops prevent basket from coming out of cart or cabinet. Baskets can be organized using adjustable dividers; order separately. Holds up to 40 pounds. Uses 4 FlexCells.

Notes

- Specify SBD8 for short adjustable divider
- Cell use with angled tilt: 6 FlexCells on 19.25"d cabinet and 8 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 23.5"w x 15.75"d x 8"h
- Interior: 21.5"w x 13.75"d x 7.75"h

■ ROAM ■ EVOLVE

Label Holder (SLH35)



Clips onto the front of all trays and baskets in a Roam cart or Evolve cabinet. Holds 2" x 3.5" cards. Comes in packages of 10.

2"H Tray (SCT2, ST2)

TEMPO ROAM EVOLVE



For use in a Roam cart or Evolve cabinet with FlexCell interior or a 3"-high Tempo drawer. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 1 FlexCell.

- Specify SCT2 for clear tray
- Specify ST2 for grey tray

Notes

- Specify STD2 for grey Preconfigured Tray with 3 Short and 2 Long Dividers
- Specify SCSD2 clear Short Tray Divider
- Specify SCLD2 clear Long Tray Divider
- Specify SSD2 grey Short Tray Divider
- Specify SLD2 grey Long Tray Divider
- 2 safety stops included in trays for Roam carts and Evolve cabinets; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 4 FlexCells on 19.25"d cabinet and 6 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 23.5"w x 15.75"d x 2"h
- Interior: 21.5"w x 13.75"d x 1.75"h

4"H Tray, Grey (ST4)

TEMPO ROAM EVOLVE



For use in a Roam cart or Evolve cabinet with FlexCell interior or a 6"-high Tempo drawer. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 2 FlexCells.

Notes

- Specify STD4 for Preconfigured Tray with 3 Short and 2 Long Dividers
- Specify SSD4 for Short Tray Divider
- Specify SLD4 for Long Tray Divider
- 2 safety stops included in trays for Roam carts and Evolve cabinets; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 5 FlexCells on 19.25"d cabinet and 7 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 23.5"w x 15.75"d x 4"h
- Interior: 21.5"w x 13.75"d x 3.75"h

8"H Tray, Grey (ST8)

TEMPO ROAM EVOLVE



For use in a Roam cart or Evolve cabinet with FlexCell interior or a 9"-high Tempo drawer. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 4 FlexCells.

Notes

- Specify STD8 for Preconfigured Tray with 1 Short and 1 Long Divider
- Specify SSD8 for Short Tray Divider
- Specify SLD8 for Long Tray Divider
- 2 safety stops included in trays for Roam carts and Evolve casework; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 6 FlexCells on 19.25"d cabinet and 8 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 23.5"w x 15.75"d x 8"h
- Interior: 21.5"w x 13.75"d x 7.75"h

2"H Narrow Tray, Grey (SNT2)

TEMPO EVOLVE



For use in Narrow 3"-high Tempo drawer or Evolve FlexCell cabinet with offset column. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 1 FlexCell.

Notes

- Specify SNTD2 for Preconfigured Tray with 2 Short and 1 Long Divider
- Specify SSND2 for clear Short Tray Divider
- Specify SSD2 for Long Tray Divider
- 2 safety stops included in trays for Evolve cabinets; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 4 FlexCells on 19.25"d cabinet and 6 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 11.75"w x 15.75"d x 2"h
- Interior: 9.75"w x 13.75"d x 1.75"h

4"H Narrow Tray, Grey (SNT4)

TEMPO EVOLVE



For use in Narrow 6"-high Tempo drawer or Evolve FlexCell cabinet with offset column. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 2 FlexCells.

Notes

- Specify SNTD4 for Preconfigured Tray with 2 Short and 1 Long Divider
- Specify SSND4 for Short Tray Divider
- Specify SSD4 for Long Tray Divider
- 2 safety stops included in trays for Evolve cabinets; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 5 FlexCells on 19.25"d cabinet and 7 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 11.75"w x 15.75"d x 4"h
- Interior: 9.75"w x 13.75"d x 3.75"h

8"H Narrow Tray, Grey (SNT8)

TEMPO EVOLVE



For use in Narrow 9"-high Tempo drawer or Evolve FlexCell cabinet with offset column. Seamless construction for easy cleaning. Trays can be organized using adjustable tray dividers; order separately. Holds up to 40 pounds. Uses 4 FlexCells.

Notes


- Specify SNTD8 for Preconfigured Tray with 1 Short and 1 Long Divider
- Specify SSND8 for Short Tray Divider
- Specify SSD8 for Long Tray Divider
- 2 safety stops included in trays for Evolve cabinets; 2 corner blocks included in trays for Tempo carts
- Cell use with angled tilt: 6 FlexCells on 19.25"d cabinet and 8 FlexCells on 27.25"d cart or cabinet

Dimensions

- Exterior: 11.75"w x 15.75"d x 8"h
- Interior: 9.75"w x 13.75"d x 7.75"h

Divider Kit (SDK3, SDK6)

■ TEMPO



For use in Tempo drawers. Includes divider guides. Additional short dividers available; order separately. Assembly required; instructions included.

- Specify SDK3 for 3"h drawer, includes 3 long and 3 short adjustable dividers
- Specify SDK6 for 6"h drawer, includes 3 long and 3 short adjustable dividers

Notes

Short Divider Package includes 5 dividers

- Specify SADP3 for Short Divider Package for 3"h drawer, dividers are 2.20"h
- Specify SADP6 for Short Divider Package for 6"h drawer, dividers are 5"h

Divider Kit for Narrow Drawer (SNDK3, SNDK6)

■ TEMPO



For use in Narrow Tempo drawers. Includes divider guides. Additional short dividers available; order separately. Assembly required; instructions included.

- Specify SNDK3 for 3"h narrow cart drawer, includes 2 long and 2 short adjustable dividers
- Specify SNDK6 for 6"h narrow cart drawer, includes 2 long and 2 short adjustable dividers

Notes

Short Divider Package, includes 5 dividers

- Specify SADP3 for Short Divider Package for 3"h drawer, dividers are 2.20"h
- Specify SADP6 for Short Divider Package for 6"h drawer, dividers are 5"h

Basic Metal Shelf (S19S1, S27S1)

■ ROAM ■ EVOLVE



Fits anywhere within a single compartment in a Roam cart or Evolve cabinet with FlexCell interior. Holds up to 25 pounds.

- Specify S19S1 for use with 19.25"d cabinet
- Specify S27S1 for use with 27.25"d cart or cabinet

Dimensions

- 23.5"w x 16"d (S19S1)
- 15.75"w x 23.75"d (S27S1)

Narrow Basic Metal Shelf (S19NS1)

■ EVOLVE



Fits in the narrow compartment of an Evolve cabinet with offset column. Holds up to 25 pounds.

Dimensions

- 11.5"w x 16"d

Metal Shelf for Cath Managers (S19M1, S19M2, S27M1, S27M2)

■ ROAM ■ EVOLVE



Fits anywhere in a Roam cart or Evolve cabinet to hold cath managers. Holds up to 50 pounds.


- Specify S19M1 for 1 compartment, 19.25"d cabinet
- Specify S19M2 for 2 compartment, 19.25"d cabinet
- Specify S27M1 for 1 compartment, 27.25"d cart or cabinet
- Specify S27M2 for 2 compartment, 27.25"d cart or cabinet

Dimensions

- 23.5"w x 16"d (S19M1)
- 32.75"w x 16"d (S19M2)
- 15.75"w x 23.75"d (S27M1)
- 32.75"w x 23.75"d (S27M2)

Adjustable Shelf Divider (S19DI, S27DI)

■ ROAM ■ EVOLVE



Adjustable and removable metal divider fits anywhere on a shelf. It is 5" high. Uses 2 FlexCells.

- Specify S19DI for use with 19.25"d shelf
- Specify S27DI for use with 27.25"d shelf

Drawers (S19D1, S19D2, S27D1, S27D2)

■ ROAM ■ EVOLVE



For use with a Roam cart or Evolve cabinet with FlexCell interior to store and organize supplies. Includes keyed lock. Holds up to 40 pounds. Uses 3 FlexCells.

- Specify S19D1 for 22.5"w x 6"h Drawer for 19.25"d, 1 compartment cabinet
- Specify S19D2 for 31.75"w x 6"h Drawer for 19.25"d, 2 compartment cabinet
- Specify S27D1 for 14.5"w x 6"h Drawer for 27.25"d, 1 compartment cart or cabinet
- Specify S27D2 for 31.75"w x 6"h Drawer for 27.25"d, 2 compartment cart or cabinet

Interior Dimensions

- 18.76"w x 15.65"d x 6.19"h (S19D1)
- 30.52"w x 15.65"d x 6.19"h (S19D2)
- 13.32"w x 23.65"d x 6.19"h (S27D1)
- 30.55"w x 23.65"d x 6.19"h (S27D2)

Cubby (ST3C, ST6C, ST9C, STN3C, STN6C, STN9C)

■ TEMPO



6"h Cubby and 9"h Cubby on Tempo cart shown

Fits in the interior of a Tempo procedure cart in place of a drawer to provide open storage organization. Comes in 3", 6", and 9" heights.

- Specify ST3C for 3"h Cubby
- Specify ST6C for 6"h Cubby
- Specify ST9C for 9"h Cubby
- Specify STN3C for 3"h Narrow Cubby
- Specify STN6C for 6" Narrow Cubby
- Specify STN9C for 9" Narrow Cubby

Note

- 9"h Cubby can only be located at bottom of cart

Interior Dimensions

- 22.75"w x 16.75"d x 2.5"h (3"h Cubby)
- 22.75"w x 16.75"d x 5.5"h (6"h Cubby)
- 22.75"w x 16.75"d x 9.5"h (9"h Cubby)

Boxed Catheter Insert (S19BC1, S19BC2, S27BC1, S27BC2, S27BC2-AT)

■ ROAM ■ EVOLVE



Boxed Catheter Insert with anti-tip mechanism shown

For use in a Roam cart or Evolve cabinet to hold catheters. Holds up to 50 pounds. Attachment hardware included. Uses 17 FlexCells.

- Specify S19BC1 for use in 19.25"d, 1 compartment cabinet. Includes 3 adjustable dividers
- Specify S19BC2 for use in 19.25"d, 2 compartment cabinet. Includes 5 adjustable dividers
- Specify S27BC1 for use in 27.25"d, 1 compartment cart or cabinet. Includes 2 adjustable dividers
- Specify S27BC2 for use in 27.25"d, 2 compartment cart or cabinet. Includes 5 adjustable dividers
- Specify S27BC2-AT for use in 27.25"d, 2 compartment cart only. Includes anti-tip mechanism and 4 adjustable dividers (2 per section)

Dimensions (Interior dimensions)

- 19BC1: 22"w x 16"d x 30-42.25"h (20.48"w x 14.48"d x 26.24-38.24"h)
- S19BC2: 31"w x 16"d x 30-42.25"h (29.48"w x 14.48"d x 26.24-38.24"h)
- S27BC1: 14.25"w x 23.50"d x 30-42.25"h (12.73"w x 21.98"d x 26.24-38.24"h)
- S27BC2: 31"w x 23.50"d x 30-42.25"h (29.48"w x 21.98"d x 26.24-38.24"h)
- S27BC2-AT: 31"w x 21.50"d x 30-42.25"h
(Each compartment - 15.15"w x 21.40"d x 26.24-38.24"h)


LED Lights (LED)

■ ROAM ■ EVOLVE

Attach to channels along the front of all Roam carts and Evolve cabinets with hinged doors to provide interior lighting. Lights are Cool White (6000k) and have an average life of 50,000 hours. Not available on Evolve upper cabinets with sloped tops.

Keyed Lock (SLK)

■ ROVER ■ TEMPO ■ ROAM ■ EVOLVE



For keyed access to a Roam cart, Tempo cart, Rover workstation, or Evolve cabinet. Includes 2 keys.

Note

- Keyed lock standard with Tempo cart and Rover workstation

Keypad Lock (SLE)

■ ROVER ■ TEMPO ■ ROAM ■ EVOLVE



For keypad access to a Roam cart, Tempo cart, Rover workstation, or Evolve cabinet.

Proximity Lock (SLP)

■ ROVER ■ TEMPO ■ ROAM ■ EVOLVE



For keycard access to a Roam cart, Tempo cart, Rover workstation, or Evolve cabinet for up to 2,000 users and administrators. Includes an LED battery indicator that notes low battery.

Dual Credential Proximity Lock (SLD)

■ ROVER ■ TEMPO ■ ROAM ■ EVOLVE



For keycard and keypad access to a Roam cart, Tempo cart, Rover workstation, or Evolve cabinet for up to 2,000 users and administrators. Includes an LED battery indicator that notes low battery.

Countertop, Compression Board with Laminate (SEC19CT, SEC27CT)

■ EVOLVE



The laminated 1.125"-thick surface fits on an Evolve base cabinet. It has a 4"-high backsplash, finished ends, and Wilsonart Frosty White (1573) laminate surface.

- Specify SEC19CT for 19.25"d cabinet; countertop is 19.75"d
- Specify SEC27CT for 27.25"d cabinet; countertop is 27.75"d
- Contact Customer Service to specify alternative laminate

Countertop, Stainless Steel (SES19CT, SES27CT)

■ EVOLVE

The stainless steel 1.5"-thick surface fits on an Evolve base cabinet. It has a 4"-high integrated backsplash.

- Specify SES19CT for 19.25"d cabinet; countertop is 19.75"d
- Specify SES27CT for 27.25"d cabinet; countertop is 27.75"d

Countertop Support Bracket (SECWSB, SESWSB)

■ EVOLVE

The bracket provides additional support for an Evolve countertop that extends past the base cabinet.

- Specify SECWSB for metal bracket used with laminate countertops
- Specify SESWSB for stainless steel bracket used with stainless steel countertops

Finished End (TFE, UFE, BFE)

■ EVOLVE



Finishes the sides of an Evolve cabinet to match the cabinet's laminate.

- Specify TFE for infill panel used with tall cabinet
- Specify UFE for infill panel used with upper cabinet
- Specify BFE for infill panel used with base cabinet

Infill Panel (TIF, UIF, BIF)

■ EVOLVE

Panel fills gaps between adjacent Evolve cabinets in a run. Panel height matches the height of the adjacent cabinet, and panel finish matches cabinets' laminate surface or brushed aluminum finish. Maximum length is 96".

- Specify TIF for infill panel used with tall cabinet
- Specify UIF for infill panel used with upper cabinet
- Specify BIF for infill panel used with base cabinet

InnerSpace

PORTFOLIO OF PRODUCTS

Solaire Medical specializes in storage solutions for all areas of your healthcare facilities. We offer the InnerSpace® portfolio of products, a cohesive system of carts, cabinets, and modular casework that improves staff and facility efficiency and an organization's bottom line.

When you purchase our products, you're buying products made in the U.S.A. and from a direct manufacturer that delivers on-time and error-free.

Contact our service team with your requests or to connect with a sales representative. We're here to help you specify the best storage products for your needs, whether a standard option or customized solution. We offer complimentary storage analysis, design services, CAD support, and logistics and installation assistance.

Call to schedule a free storage consultation.

solaire
MEDICAL

Solaire Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501
888-435-2256 phone | 616-743-5944 fax | solaiomedical.com



Roam
Supply Carts



Rover
Workstations



Tempo
Procedure Carts



Evolve
Cabinets



Evolve
Casework