## ECP Cart with Drawers and 2 Center Columns (SRSDCECP)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 3 ECP (extending cardiovascular procedure) Cart with Drawers and 2 Center Columns has a stainless steel top that extends both left and right to increase surface space, 2 full-width drawers, and 5 FlexCells in each of the cart's 3 compartments. The drawers have full-length extension glides and integral pulls. The cart comes in a Soft White powder-coated paint finish.







## **Materials and Finishes**

Base structure: **Powder-coated steel**Top structure: **304 Stainless Steel** 

Corner extrusions: Clear anodized aluminum

Corner bumpers: Thermoplastic polymer (ABS plastic)
FlexCell panels: Thermoplastic polymer (ABS plastic)

Drawers: **Powder-coated steel**Casters: **Nylon and steel** 

**Dimensions** (nominal) and Weight

Exterior: 57"w x 28.75"d x 36"h

Interior: 14.5" w (each compartment) x 23.75" d x 12.75"h

Weight: 330 lbs

## **Interior Configuration**

- 15 FlexCells, 5 per compartment
- 2 full-width drawers

## **Options**

Interior and side-mounted exterior accessories

