
## ECP Cart with Drawers and Dividers (SRSDECP)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 3 ECP (extending cardiovascular procedure) Cart with Drawers and Dividers has a stainless steel top that extends both left and right to increase surface space, 2 full-width drawers, and 1 shelf with 5 adjustable dividers. The drawers have full-length extension glides and integral pulls. The cart comes in a Soft White powder-coated paint finish.





Storage solutions for healthcare





Base structure: **Powder-coated steel** Top structure: **304 Stainless Steel** Corner extrusions: **Clear anodized aluminum** Corner bumpers: **Thermoplastic polymer (ABS plastic)** Shelf: **Powder-coated steel** Drawers: **Powder-coated steel** Casters: **Nylon and steel** 

Dimensions (nominal) and Weight

Exterior: 57"w x 28.75"d x 36"h Interior: 53"w x 26"d x 12.75"h Weight: 318 lbs

**Interior Configuration** 

· Shelf with 5 adjustable dividers and 2 full-width drawers

**Options** 

Interior and side-mounted exterior accessories



Storage solutions for healthcare

Solaire Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501 888-435-2256 phone | 616-743-5944 fax | solairemedical.com