


Evolve TEE Probe Cabinet with Roll-Top Door

Evolve modular cabinets and casework create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts.

The 19.25"-deep, 92"-high TEE probe cabinet comes in 2 widths. It has a counterbalanced roll-top door, is vented on the top and bottom, and includes a removable drip tray. Material choices include compression board with a melamine or specified laminate surface or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. Cabinet has a 4"-high toe kick.



Evolve TEE Probe Cabinet with Roll-Top Door



Dimensions (Dimensions are nominal)

Exterior (interior):

26.75"w x 19.25"d x 92"h (25.25"w x 17"d x 79.5"h)

36"w x 19.25"d x 92"h (34.5"w x 17"d x 79.5"h) (shown)

92"h cabinet with 6"h toe kick is 94"h

Interior Configuration:

26.75"w cabinet holds 3 TEE probes

36"w cabinet holds 5 TEE probes

Options

- Optional Fan Ventilation Unit available
- Keyed or keyless lock options available
- Optional sloped top available
- Optional 6"h toe kick available

Product Number Reference



Product Shown:
SEC1936PRTT

S E C 1 9 3 6 P R T T

Product Line

Material

Depth

Width

Cabinet
Type

Door Type

Hinge
Location

Cabinet Configuration



Product Numbers

26.75"w x 19.25"d x 92"h

Compression board: **SEC1927PRTT** (shown)

AireCore: SEA1927PRTT

AireCore and Brushed Aluminum: SEB1927PRTT

Stainless steel: SES1927PRTT



36"w x 19.25"d x 92"h

Compression board: **SEC1936PRTT** (shown)

AireCore: SEA1936PRTT

AireCore and Brushed Aluminum: SEB1936PRTT

Stainless steel: SES1936PRTT

Evolve TEE Probe Cabinet with Roll-Top Door

Materials and Finishes

Cabinet shell:

Compression board: Thermally fused laminate over 0.75"-thick compression board; melamine surface or specified Wilsonart, Formica, or Pionite standard laminate

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish

Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

Edge banding:

Compression board and AireCore: White ABS

AireCore with Brushed Aluminum finish: Brushed Aluminum acrylic

FlexCell panel:

Injection-molded ABS

Roll-top door and recessed channel:

Polypropylene



Storage solutions for healthcare

Solairé Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501
888-435-2256 phone | 616-743-5944 fax | solairomedical.com

© 2018 Solairé Medical, Grand Rapids, MI 49501
© InnerSpace is a registered trademark of Solairé Medical, Grand Rapids, MI 49501