

**PH-31**

## MRI Quality Assurance Phantom MHR

*This QA phantom for MRI allows to evaluate the slice thickness, spatial resolution, uniformity and geometric distortion as well as contrast. Complies with NEMA standards.*



### Specifications


Set Includes:  
1 phantom unit A  
1 phantom unit B  
1 set of nickel chloride solution  
1 storage case

Size:  
22 dia. x 14cm/8.7 dia. x 5.5 in

**PH-26**

## ORINS Thyroid Phantom ITS

*This phantom is Oak Ridge Institute for Nuclear Studies (ORINS). For measurement of RI thyroid uptake.*



### Specifications

Set includes:  
1 storage case

Size:  
12.5 (dia) x 12.5 (H) cm/4.9 x 4.9 (H) in.

**PH-32**

## MRI Quality Assurance Phantom JMR

*SN ratio, uniformity, slice thickness, spatial resolution, geometric distortion can be evaluated.*



**Conforming to JIS Z 4924**

### Specifications

Set Includes:  
1 phantom unit A  
1 phantom unit B  
1 set of nickel chloride solution  
1 storage case

Size:  
22 dia. x 14cm/8.7 dia. x 5.5 in

**PH-32B** NEW

## MRI Quality Assurance Phantom JMR II

*This phantom was developed corresponding to JIS Z 4924 revised in 2016.*

*Both measurement of the ghost and evaluation of the image contrast are added to test items. The size of this phantom is modified to the human head size.*



### Specifications

Set includes;  
1 phantom  
1 storage case

Size  
Size: 180 (dia) X 130(H)