


Roam 2 Long Catheter Cart (SR2GC)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 2 long catheter cart comes with 7 9-hook cath managers. Cath managers are permanently mounted to the top of the cart and can be adjusted horizontally. The cart has a lightweight aluminum case, durable bumper, easy-grip bars for maneuvering the cart, 5" medical-grade casters, and comes in a Soft White powder-coated paint finish.



Roam 2 Long Catheter Cart (SR2GC)



Materials and Finishes

Base and top structure: **Powder-coated steel**
Corner extrusions: **Clear anodized aluminum**
Corner bumpers and corner trim: **Thermoplastic polymer (ABS plastic)**
Hinges: **Steel**
Door pulls: **Dull chrome-plated steel**
Push handles: **Stainless steel**
Glass doors: **Clear anodized aluminum and tempered glass**
Casters: **Nylon and steel**

Dimensions (nominal)

Exterior: **39"w x 28.75"d x 80.5"h**
Interior: **34.75"w x 23.75"d x 68.5"h**
Weight: **215 lbs**
Interior Configuration: **7 9-hook cath managers; 64.8" hanging length**

Options

- Keyless or keyed lock
- 5" sloped top
- Five paint colors and range of laminates

