


Roam 2 Scope Cart (SR2GS)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 2 scope cart holds 8 endoscopes, up to 60" in length. It has tempered glass doors, a cushioned back panel to protect scope optics, a ventilated bottom, and a removable drip tray. Scope managers and cord retainers are adjustable and keep contents organized and elevated off the cart floor. The cart comes in a Soft White powder-coated paint finish.



Roam 2 Scope Cart (SR2GS)



Materials and Finishes

Base and top structure: **Powder-coated steel**

Corner extrusions: **Clear anodized aluminum**

Corner bumpers and corner trim: **Thermoplastic polymer (ABS plastic)**

Hinges: **Steel**

Door pulls: **Dull chrome-plated steel**

Push handles: **Stainless steel**

Glass doors: **Clear anodized aluminum and tempered glass**

Casters: **Nylon and steel**

Dimensions (nominal)

Exterior: **39"w x 28.75"d x 75.25"h**

Interior: **34.75"w x 23.75"d x 64.25"h**

Weight: **200 lbs**

Interior Configuration: **Holds up to 8 scopes, maximum length: 60"**

Options

- HEPA fan unit
- Keyless or keyed lock
- 5" sloped top
- Five paint colors and range of laminates

