


Roam 3 Supply Cart with Roll-Top Doors (SR3R)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 3 cart comes with two counterbalanced roll-top doors, three interior compartments, and in a Soft White powder-coated paint finish.




Roam 3 Supply Cart with Roll-Top Doors (SR3R)



Materials and Finishes

Base and top structure: **Powder-coated steel**
Corner extrusions: **Clear anodized aluminum**
Roll-top track extrusion: **Clear anodized aluminum**
Corner bumpers and corner trim: **Thermoplastic polymer (ABS plastic)**
FlexCell panels: **Thermoplastic polymer (ABS plastic)**
Push handles: **Stainless steel**
Roll-top door slats and recessed channel: **Polypropylene**
Casters: **Nylon and steel**



Compartment configurations

Dimensions (nominal)

Exterior: **58.75"w x 28.75"d x 81"h**
Interior: **14.5"w x 23.75"d x 64.25"h** (single compartment)
31.75"w x 23.75"d x 64.25"h (double compartment)
Weight: **288 lbs**
Interior Configuration: **75 FlexCells, 25 per compartment**

Options

- With single compartment on the left and two center columns, with right center column, or with left center column
- Electronic or keyed lock
- Interior and side-mounted exterior accessories
- Five paint colors and range of laminates

