



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

9941 ULTRALITE

MODEL: [RE-9941/9941A/994153]

DESCRIPTION:

The 9941 Ultralites come standard with a rubberized nose bridge. They are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields. The 9941 is a popular frame and ideal for those who have a high power prescription. Available in 2 sizes and optional adjustable nose bridge for regular sized frames (9941A).

DESIGN:

Model: 9941 (regular) / 9941A (nose bridge)
Frame style: Wraparound
Front lens: 0.75mm Pb Glass
Frame Material: Plastic
Weight: 64g
Frame size: 55-22-130
Side Shields: 0.50mm Pb strips (optional)

Model: 994153 (small)
Frame style: Wraparound
Front lens: 0.75mm Pb Glass
Frame Material: Plastic
Weight: 60g
Frame size: 53-22-130
Side Shields: 0.50mm Pb strips (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00
Minimum Progressive PD: 64
Base Curve: 6-4

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00
Minimum Progressive PD: 58
Base Curve: 6-4

FEATURES:

Adjustable Nose Pads (9941A), Rubberized Temples

COLORS:

Silver (SIL), Black (BLK), Blue (BLU), Red (RED)



SIL



BLK



BLU



RED