

Ceiling column mounting plate

300/P/004 (Installation aid for securing a Kenex ceiling column)

An effective installation aid

This mounting plate can be fixed below Unistrut® or equivalent ceiling channels, which form part of a suspended ceiling. It provides a simple and effective way of installing a ceiling column. **Note:** The load bearing capacity of the ceiling channels must be checked by structural analysis and confirmed by means of an acceptable certificate.

The plate measures: 735 mm x 340 mm x 38 mm. The installer must ensure that the two ceiling channels are at 650 mm +/- 45 mm apart, and are horizontal and parallel.



Weights and dimensions

Size: 735 mm x 340 mm x 38 mm

Weight: 15 kg

Part number: 300/P/004