

Featherweight® Model FWSTII Single Tower Endoscopy Cart

Innovation Designed Around You

**PHS
WEST
INC.**

Technical Specifications

If your requirement for a small footprint necessitates a Single Short but the available area for your equipment and supplies do not fit, the Featherweight® Single Tower Endoscopy Cart give you the vertical space you need.



Ergo-Express® Model FWSTII | Single Tower Endoscopy Carts

Specifications

- Model Number: FWSTII
- Dimensions:
 - Depth: 24"
 - Width: 25"
 - Height: 5=57.5" | 8=60.5"
 - Weight: Dependent upon cart size
 - Shipping Weight: Dependent upon cart size
- Extrusion Dimensions and Materials
 - 30mm x 60mm Aluminum
- Caster Wheels
 - Plated Carbon Steel
- Bumpers
 - Vinyl
- Drawer Options
 - Sizes: 4" | 6" | 8" | 12"
 - Width: 20.25"
 - Height: 3.5" | 5.25" | 7.125" | 10.75"
 - Depth: 18.375"
- Materials
 - Plastic Panels
 - ABS
 - Shelves and Deck
 - Aluminum
 - Work Surface
 - PVC
- Three Year Limited Warranty
- [24/7-365 Technical and Operational Support for the life of the product](#)

Additional technical and operational specifications are available for custom designs.

