




Mobile DR detector and cassette holder

1305/8/1 (Suitable for portable DR Panels and CR or film cassettes)



Features

- > Ideal for imaging examinations i.e. mammography.
- > Unlimited orientation for universal positioning possibilities.
- > Telescopic width adjustable from 24.5 cm to 53.5 cm.
- > Holds full leg/spine image receptors of any length.
- > Not counterbalanced, but with easy height adjustment from floor level to 142 cm.
- > Overhangs table or trolley by up to 49 cm.
- > Accepts image receptor/cassette up to 30 mm thick of any length.
- > Heavy base and low centre of gravity gives added stability.
- > Weight 52.5 kg.