




Patients imaging chair

1501/D (an ideal chair for a range of imaging procedures)



Features

- > Ideal for imaging examinations i.e. mammography.
- > Hydraulic operation alters seat height between 60 cm & 75 cm.
- > Height adjustment made with ease regardless of the weight of the seated patient. Maximum weight 128 kg (282 lbs).
- > A single lever can raise, lower or lock the rotation of the seat
- > The back and arm rests are held securely yet can be removed & replaced in seconds.
- > Each castor has a rolling and swivel brake, so any locked pair will prevent movement.

Weights, dimensions & ordering information

Model: **1501/D**

Weight: **24 kg**

