

Wall rack

312/W-TOP

Wall rack for upper removable top shields

This wall rack provides a secure storage location for the removable upper tops, used on Kenex lower body radiation shields. Maximum configuration: one 57 cm wide top and two 27 cm wide tops, or four 27 cm tops.





Shields not included





Technical information



Any top shield used in conjunction with the Kenex range of lower body radiation shields (types 311 & 312), easily drops into slots machined into this strong, plastic wall rack. The compact size enables the tops to be located within a small area on a wall.

Weights & dimensions

Width: 600 mm Depth: 146 mm Height: 108 mm Weight: 1 kg

Supplementary products


Wall rack for lower body radiation shield: 312/W-LB Effective storage for any size Kenex lower body radiation shield.





Mobile rack for table shielding: 1342/A2

Effective mobile storage combinations for lower body table shields (types 311 & 312), removable Kenex top shields and optional head end shields (types 313).



Please note shielding is not included and is shown here for display purposes only. Further information on these products is available on our website and by downloadable data sheet.

