

# EEBP

Optimize  
E-Prime experiment  
communication with  
Brain Products  
EEG systems



## Features:


- Brain Products toolbox of package calls in E-Studio
- Drag and drop package calls give the ability to:
  - Start, stop, pause and continue your EEG recording
  - Send markers at specified events (e.g., onset of stimulus or response)
  - Send customized markers based on stimuli type for simpler data analysis
  - Send marker of response accuracy
  - Verify recording has started
  - Receive notification of EEG errors
  - Load Recorder Workspaces
  - Set Recorder Monitoring, Impedance, and Test Signal modes

## Includes:

- User manual with experiment tutorial
- Sample experiments for Trigger Box, parallel port, and Chronos communication
- One year of support provided by PST

## Requirements:

- E-Prime 3.0
- BrainVision Recorder and Remote Control Server



412.449.0078



sales@pstnet.com



412.449.0079



www.pstnet.com

