**Optimize E-Prime** experiment communication with **EGI EEG systems** 



## Features:

- Maintain time synchronization between E-Prime and Net Station
- EGI experiment control interface (ECI) available to integrate with other experiment control systems
- Net Station toolbox of package calls in E-Studio
- Drag and drop package calls give the ability to:
  - · Start and stop your EEG recording
  - · Load, start, and stop pulse trains
  - · Send markers based on onset of objects and trial condition for simpler data analysis


## Includes:

- User manual with experiment tutorial
- Samples and Tutorials
- All Package Calls needed for seamless integration
- E-Prime support and EGI support

**Requirements:** 

- E-Prime 3.0
- EGI Net Station software

























412.449.0079 ( www.pstnet.com