



MOBILE MULTIPLE APRON HANGER WITH 10 SWING ARMS MODEL RAR-M10S-HB



RAR-M10S-HB
(10 Aprons)

DESCRIPTION

The Phillips Mobile Multiple Apron Hanger with 10 Swing Arms has a generous, sturdy base and solid steel construction, making it virtually untippable. It can safely hold up to 10 lead aprons, and its Mobile Pedestal System allows it to be move between radiation departments, labs, operation rooms, and more.

FEATURES

- *Stable, Sturdy Base*
- *Holds Up to 10 Aprons*
- *Mobile Pedestal System*

SPECS

- *Overall Dimensions: 24" D x 20" W x 49.5" H*
- *Weight: 60 lbs.*
- *Arms Extend a Maximum of 22.5" in Depth*



Phillips Safety Products
123 Lincoln Blvd.
Middlesex, NJ 08846

Contact Us
1-866-575-1307
info@phillips-safety.com